PULSOSIMETRO
PORTATILE

PL101

IT EN FR ES

PULSOSIMETRO PORTATILE

IT

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application site regularly. Vary the sensor application site and check skin integrity, circulatory status and correct orientation at least every 4 hours.

- Inaccurate measurements could be caused by autoclaving, ethylene oxide sterilization, immersion of the device in liquids, significant levels of dysfunctional hemoglobins (such as carboxy-hemoglobin or methemoglobin), intravascular dyes such as indocyanine green or methylene blue.
- SpO₂ measurements could be negatively affected in the presence of high ambient light (avoid direct exposure of the sensor to sunlight), excessive patient movements, vein pulsations, placement of the oximeter on a finger and simultaneous presence of a pressure cuff on the arm, from the presence of an arterial catheter or an intravascular line, in patients with hypotension, severe vasoconstriction, with severe anemia or hypothermia, in the presence of cardiac arrest or in a state of shock.
- The presence of nail polish or false nails can cause inaccurate readings.

FEATURES

- OLED display
- Product adopts double color OLED display can show the six different display mode
- Low-power consumption, continuously work for 40 hours
- Low voltage indicator
- In the absence of signals, the product will be in after 10 seconds to enter standby state
- Visual and auditory prompt function.
- Small in volume, light in weight, and convenient to carry
- Product Operation Scope
- The fingertip oximeter can be used to measure human haemoglobin saturation and heart rate through finger. The product applies to using in family, hospital (including clinical use in internist/surgery, anaesthesia, paediatrics, intensive care, etc.), oxygen club, social medical organizations, physical care in sports (It can be used before or after sports. Operation in sport procedure is not recommended). It is also applicable to enthusiasts on mountaineering, patients (convalescents at home or those need first aid treatment), elders over 60, those work more than 12 hours, sporters and those work in the hermetic circumstance, etc. The product is not suitable to monitor patient continuously.

OPERATION INSTRUCTIONS


- Hold the lock button on the battery cover, push the battery cover down at the same time, the battery cover is opened.
- Install two AAA batteries into battery cassette before covering its cover.
- Plug one finger into rubber hole of the oximeter (it is best to plug the finger thoroughly) before releasing the clamp with the nail upwards.
- Press button on the front panel;
- Don't tremble your finger when the oximeter is working. It is recommended to stay still while using your oximeter.
- Press the button on the front panel, if we want change display direction;
- Read relevant datum from display screen.
- The instrument has the function of sleep, no signal will enter standby state of sleep;
- Please replace new batteries when OLED indicates the batteries are in low power.

When plugging your finger into the oximeter, your nail surface must be upward.

Note: Please use medical alcohol to clean the rubber before each test and clean the tested finger with alcohol before and after the test. (The rubber inside of the oximeter adopts medical rubber, which has no toxin, no harm, and brings no side effect such as allergy to the skin).

BRIEF DESCRIPTION OF FRONT PANEL

When the measurement is displayed, pressing the ON / OFF key changes the orientation of the display to make reading easier: in positions 1 and 2 the Plethysmographic Wave (pulse wave) is also visible.

DETAILED DESCRIPTIONS OF PRODUCT FUNCTIONS:

- Display type: OLED
- SpO₂ measurement range: 70% -99%
Precision: for SpO₂ values between 70% - 99% ± 2%; for values (<70%) not specified Resolution: ± 1%
- Measurement range PR: 30BPM-240BPM
Accuracy: ± 1BPM or ± 1% (greater)

4. Parameters of LED probe

	wavelength	radiation power
RED	660±2nm	1.8mW
IR	905±10nm	2.0mW

This parameter can be especially useful to clinicians.

- Power supply: two 1.5V alkaline batteries (AAA)
- Supply current: <30mA
- Automatic shutdown: after 10 seconds without signal.
- Environmental characteristics:

Working temperature: 5°C-40°C
Storage temperature: -10°C-40°C

Working humidity: 15%-80%

Environmental conservation humidity: 10%-80%

Working pressure: 70kPa -106kPa

The device complies with IEC 60601-1-2 for EM emissions. The materials which the user can come into contact with have no toxicity and no action on the tissues and are in accordance with ISO 10931-1, ISO 10935-3 and ISO 10933-10.

Subject to technical modifications.

CLASSIFICATION

Management Class for Medical Devices: equipment
Anti-electric Shock Type: Internally powered equipment
Anti-electric Shock Degree: Type BF equipment

MAINTENANCE AND PRESERVATION

- Replace the batteries timely when low voltage lamp is on.
- Clean the surface of fingertip oximeter before it is used to diagnose patients.
- Remove the batteries inside if you will not operate the oximeter for a long time.
- It would be better to preserve the product in -10~40 (14-104) and humidity is 10%-80%.
- It is recommended that the product should be kept dry anytime. A wet ambient might affect its lifetime and even damage the product.
- Please follow the law of the local government to deal with used batteries; subjected to WEEE Directive, dispose of according to local laws.

PRODUCT ACCESSORIES

- Lanyard and bag for carrying
- Two 1.5V AAA alkaline batteries
- User Manual

GUIDANCE AND MANUFACTURE'S DECLARATION - ELECTROMAGNETIC RADIATION - FOR OTHER EQUIPMENTS AND SYSTEMS

The Pulse Oximeter is designed to be used in specified electromagnetic environment. Users of the Pulse Oximeter must use it in the following environments.

Radiation Test	Compliance	Electromagnetic environment-guidance
RF interference	Group 1	RF signal of Pulse Oximeter is simply created by its internal function. Therefore, its RF interference is very low and is not likely to cause any interference to nearby electronic equipment.

POSSIBLE PROBLEMS AND RESOLUTIONS

Problem	Possible reason	Solution
SpO ₂ or PR can not be shown normally	1. Finger is not plugged correctly 2. Patient's Oxygen-hemoglobin value is too low to be measured	1. Retry by plugging the finger 2. Try more times. If you can make sure there is no problem in the product, please go to hospital timely for exact diagnosis
SpO ₂ or PR is shown unsteadily	1. The finger might not be plugged deep enough 2. Finger is trembling or the patient is on movement status	1. Retry by plugging the finger 2. Please remain at rest
The oximeter cannot be turned on	1. Inadequate power or power off 2. Batteries might be installed incorrectly 3. The ox	

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